

# ENdynamics Ltd

Enterprise Centre Washington St. Bolton BL3 5EY Tel: 01204 537792 Email: sales@endynamics.co.uk

## Manual Wheelchair – The Eclipse

Nothing less than the world's best Bariatric wheelchair! Its height, depth and width can be adjusted. In short it follows the users weight gain or weight loss, or can be re-built for another user.

The unique wheelchair frame design ensures that the greater part of the weight is taken up on the rear wheels. This lightens the steering and reduces energy loss.

Its seat depth and seat to floor height are adjustable, seat and back can be angled, allowing greater variation in adapting to the individual user.

The heavy-duty steel frame has a profile that ensures maximum long-lasting strength and rigidity.

With its removable back, arm, footrests and quick-release wheels the Eclipse can be dismantled to fit into smaller cars quickly and easily.

Built for large individuals, this wheelchair will cope with the extra weight strain with style and ease.

- Easy to drive, easy to push
- Specially designed for bariatric users
- Anti-tip in front and at back



## Accessories

- Strapback
- Canestyle pushbar
- Reclining back
- IV holder
- O2 tank holder
- Head support
- Elevation footrest
- Full-size footplate flip-up

Maximum weight-bearing capacity	270/450 kg
Seat width	56 - 115 cm
Seat depths	41 - 56 cm
Seat to floor height	36 - 51 cm
Total width	18 cm + seat width
Height of arm rest	29 - 34 cm
Weight	27 kg